

# AccuPen™ Dispensing Device

Measured Dose for Measured Results

Our pharmacy is proud to introduce you to the patient-friendly AccuPen, an accurate dispensing device for your topical compounded preparations. The AccuPen is engineered to be an airless metered dose pump, with a unique design that allows you to conveniently transport your medication. Unique to the AccuPen is a clear window that shows the amount of preparation remaining in the pen.

## Highlights and Key Characteristics

**Unique Locking Feature!** Turn the End Cap clockwise to 'lock' the device – ideal for traveling, carrying in bags, purses, etc.!




**Simple to Use:** First, remove the cap. To use the AccuPen, simply turn the Push-Button End Cap counter-clockwise to unlock the actuator, and then push the Push-Button End Cap in the same manner you would push an ink pen. A metered volume of the topical preparation will then be pumped out of the tip! It's that easy!

**How It's Made:** Sturdy design featuring a polypropylene (PP) body, a low density polyethylene (LDPE) diaphragm, and a PP, LDPE, stainless steel pump.

Ask your pharmacist for the AccuPen for easy and convenient dispensing of your topical medications.

Have questions? Ask the pharmacy staff – we're here to help!



distributed by  
 PCCA